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
Instructions for use for Dorsolumbar Spinal Brace (Taylor)

CODES: 10.06.020

EN

Instructions for use

2S Dorsolumbar Spinal Brace (Taylor)

SPECIFICATIONS / INTENDED USE

If you are diagnosed with a spinal problem or have had spinal surgery your physician may prescribe the use of an external support structure, such as the Dorsolumbar Brace. Such spinal problems may include spinal deformity, vertebral fractures or collapse, osteoporosis of the spine, congenital abnormalities of the lumbar sacral region and spondylolisthesis. Some spinal surgeries where the use of a brace may be prescribed are a spinal fusion, laminectomy or discectomy. The brace is designed to prevent excessive bending, control the upper and lower regions of the spinal column and support your spine whilst standing, sitting and walking. It also offers a safe, non invasive way to prevent future problems or to help you heal from a current condition.

APPLICATION INSTRUCTIONS

MEASUREMENT & SIZING

Proper sizing is important. Measure around your low back region to determine the correct size. If you have a large mid section, consider getting the next size up so that you will be comfortable.

S	M	L	XL	XXL
70-80cm	90-90cm	90-100cm	100-110cm	110-120cm

XXXL	XXXXL
120-135cm	135-150cm

APPLICATION

The detailed application instructions of the spinal brace will be provided by your physician. These instructions need to include the number of hours per day, the estimated time period for wearing the device (weeks, months, lifetime), how tight the device needs to be fitted, and whether or not short periods of "relaxed" fitting pressure are permissible during the active wearing time.

Since the brace is designed to limit forward bending, the longer the system is in the front, the more effective it becomes. But because our bodies shorten when we sit, the brace usually needs to accommodate the sitting position. When seated and sitting upright with the spine as straight as possible, the front section should be no higher than the bottom of the breast and no lower than the pubic bone (just above the private area.) The back (posterior) part is designed to provide a counter pressure and has been located by your orthotist according to your pathology.

The brace is made of 2 separated parts that get attached together by strong Velcro strips in the back side.

- A. A strong, tall (32cm high) back supporting belt
- B. A vertical part featuring 2 long supportive metal strips, adjustable in height, attaching to the belt in the back side with strong Velcro and 2 extra long elastic bands that go over your shoulders, under your armpits and around your waist securing the two parts together and helping adjust the level of pressure and tightness you apply to your brace.


The first time you wear your brace, it is better that your physician helps you put it on and also they align and tighten the back part of the brace the way it should be according to your condition.

- Centre the garment (belt with supportive metal back strips attached) on your spine. Note: Ensure the fabric fasteners at the ends of the outer layer are not yet fastened to the inner layer
- Make sure the width of the belt covers the lumbar (low back) and sacral (tail bone) area.
- Do not suck your stomach in all the way, and do not let your midsection (belly) hang out all the way before putting on the belt. Start with a neutral position in which you are comfortable

- Wrap the belt around to the front. When you Velcro the first part of the belt, make sure the two surfaces are matched up and you have a solid foundation from overlap. You do not want to have a criss-cross pattern, or your belt will ride up and not stay in the proper position.


- Your belt has extra support flaps, use those for extra support around your waist and bring them to the front and Velcro fully.

- Pull the elastic bands over your shoulders and under your armpits. Make sure the protective pads are placed under your armpits in order to avoid discomfort by the tightened elastic band in this area.

- Cross the elastic bands in the back and bringing them to the front (you might need assistance for this step)

- Secure the elastic bands in the front of the garment with the Velcro. Pull as much as you feel comfortable and in order for your upper back to become as straight as needed.

Once all the Velcro straps are in place, check the position is comfortable over the hips and waist and that the shoulder straps in the correct position and tightened.

You should not feel restricted in taking a deep breath. If you are unable to take a deep breath then the brace may be fitted too tightly, and needs adjustment.

You can now sit and move around within the guidelines of your doctor

COMPOSITION

Elastic belts are made of high quality elastic suppression fabric. Cotton 40%, Polyamide 45%, Elastan 15%

GENERAL WARNINGS

- Wear brace as directed by your physician.
- Brace should be applied before getting out of bed unless instructed otherwise
- Weight changes will affect brace fit. Call your physician to schedule a brace check/adjustment.
- Always wear an under-shirt or similar garment to help protect the skin. This must 100% cotton, tight and without wrinkles. Change this shirt daily or more frequently if perspiration is a problem.
- Do not place anything aside from light clothing under your brace as it may cause a pressure area.
- Spinal braces are most comfortable and effective when worn snugly. This will help reduce slippage and rubbing.
- Re-tighten straps throughout the day to maintain proper snugness.
- Your physician will decide the duration of time for wearing the brace.
- Your spinal brace should not be worn in the shower. If your physician requires you to wear the brace in the shower, you will need to be fitted with a special brace for showering.
- If any of the metal bands are uncomfortable or digging in, do not try to modify the frame yourself. If you are uncomfortable, consult your physiotherapist.

Side effects harmful to health have not been reported to date. Please consult your clinician if you have the following:

- Skin disorders or injuries in the relevant part of the body - particularly if inflammation, raised scars with swelling, redness or excessive heat buildup is present.
- Notify your orthotist of any redness or soreness that does not go away during the period that the orthosis is off. If skin breaks down, remove the orthosis and do not reapply until instructed to do so by your orthotist or physician. While the orthosis is not being worn, treatment is stopped. Therefore, it is very important to notify your orthotist in a timely manner when skin problems occur.
- Impaired sensation.
- Heart or lung conditions.

The brace should only be used for its intended purpose in the manner described in this manual. No liability is accepted for improper use of the product.

- Sharing of the brace is not recommended.
- Please keep the brace out of the reach of children. It is not a toy.
- Should you notice any unusual changes (e.g. increased symptoms), please contact your healthcare specialist without delay.
- Avoid the use of ointments/ lotions or any substances including grease or acids during use of the brace.
- Check the product's composition for possible sensitivity/allergic reaction to any of the product's materials.

Daily Activities

Although the brace will support you, the level of activity you can do will be dictated by your condition. Your consultant and physiotherapist will lead you in this.

Walking

When wearing the brace you may walk normally, including use the stairs.

Sleeping

Your consultant will give you specific instructions as to whether you need to wear your brace when sleeping or if you can leave it off at night.

Showering

You can take sponge baths on your bed, after carefully removing your brace. After drying yourself you can put on your brace again, following your physicians instructions. Still it is very possible (depends on your condition) that you should wear your brace at all times, including during showering. If your physician requires you to wear the brace in the shower, you should better be fitted with a special brace for showering. If this is not possible, then you should shortly remove the brace, very carefully dry its fabric parts, and then put it on again.

Driving

It is best to avoid driving while wearing your brace. The use of safety belts and air bags can also affect your condition. You must get permission from your consultant first as there may be restrictions depending on your surgery or your overall condition.

Travel

You may travel as you wish with the brace once your condition allows and you have your consultant's permission.

WASHING / CLEANING INSTRUCTIONS

Close the Velcro bands, and wash the product with the Velcro straps hermetically shut.
Remove the metallic supportive accessories where existing.
Never machine wash. Hand wash in water at 30°C. Do not wring to remove water. This is to protect the internal components from damage.
Do not dry clean. Do not use bleach. Do not tumble-dry
Always dry lying flat. Do not dry with heat. (E.g. over a radiator). Do not iron
Make sure that the product is 100% dried before use.

IMPORTANT NOTICE!

WHILE EVERY EFFORT HAS BEEN MADE TO OBTAIN MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THIS DEVICE IS ONLY ONE ELEMENT IN THE OVERALL TREATMENT PROGRAM ADMINISTERED BY A MEDICAL PROFESSIONAL.
2S DISCLAIMS ANY WARRANTY, EXPRESSED OR IMPLIED THAT THESE PRODUCTS PREVENT INJURIES.

IN CASE THE PRODUCT HAS A DEFECT PLEASE CONTACT YOUR LOCAL SUPPLIER, STORE OR PHARMACY. IN CASE OF A MANUFACTURING MALLFUNCTION YOU MAY RETURN THE PRODUCT WITH YOUR PURCHASE DOCUMENT WITHIN 5 DAYS FROM THE DAY OF THE PURCHASE.